

All of the home's lifestyle sources should be removed and stored outside of the areas to be sampled. We are looking for TCE/PCE, these sources include but are not limited to:

- Any products with TCE or PCE
- Recently "Dry Cleaned" items
- Chlorinated solvents/degreasers
- Electronics cleaners
- Penetrating oils
- Rust buster/Liquid wrench
- Brake parts cleaner
- Carburetor cleaner
- Belt Dressing
- Wire drier
- Tire shine
- Plumber's goop adhesive/sealant
- Gun cleaners
- Recently cleaned guns
- WD 40
- Paint/spot removers
- Some metal polishes
- Carpet cleaners
- Some art & crafts adhesives
- Spray film cleaner
- Spray clear fixative
- Some fabric protectors
- Boot protectors/sealants/waterproofing
- Spray on snow
- Computer dust off
- White-out
- E6000 Crafter's Glue
- Mace

As a general rule the older the product the more likely it is to contain chlorinated compounds. Old rusted spray (non-paint) cans should be removed as well. Oddly enough the older products sometimes list the "chlorinated" contents. You can checkout the following links as well.

<http://householdproducts.nlm.nih.gov/cgi-bin/household/brands?tbl=chem&id=141&query=tetrachloroethylene&searchas=TblChemicals1&prodcats=all>

<http://householdproducts.nlm.nih.gov/cgi-bin/household/brands?tbl=chem&id=266&query=trichloroethylene&searchas=TblChemicals1&prodcats=all>